

Date: 4/17/2026 3:01 AM **Total Pallets:** **Print Sort:** Bin Location
Total Orders: 9 **Total Cartons:** 27 **Total Pick Labels:** 262
Picker: **Pick Date:** **Batch No.:** B20260417030051
Zone: **Total Item Qty:** 263 **Pick Type:** FC
Bin Category: **Conveyable:** **Schedule Ship Date:** 4/15/2026
Lift Type: **Rework:** N **Expedite:** True

Cust. Specific Attr.: **Individual:**
Consolidated: No **Smart Picking:** No
Pick List No.: B20260417030051FC001



| Item No. | Description | Rev. | Rev. Desc. | Kit Item No. | Loc. | Bin Loc. | Qty to Pick |
|------------|---|------|------------|--------------|------|----------|-------------|
| PU3033AH-S | 7.5Fr Ureteroscope, DISS, Standard Deflection | | | | GA2 | BB48-01A | 20 |
| PGA | Single use grasping forceps | | | | GA2 | BB46-02A | 20 |
| PU3033A-S | 7.5Fr Ureteroscope, Standard Deflection | | | | GA2 | BB39-02A | 10 |
| 205W11F40 | Suction Access Sheath 11/13Fr 40cm | | | | GA2 | BB38-01A | 10 |
| 104F150PS | Collection Cannister for FANS | | | | GA2 | BB37-02A | 20 |
| PC200-AR | Single Use Cystoscope - Reverse Deflection | | | | GA2 | BA29-01A | 100 |
| PU3033A-R | 7.5Fr Ureteroscope, Reverse Deflection | | | | GA2 | BA31-02A | 30 |
| AB01 | Bracket | | | | GA2 | BA44-01A | 2 |
| PU400A-S | 9.2Fr Ureteroscope, DISS, Standard Deflection | | | | GA2 | BB46-01C | 10 |
| PC200-AS | Single Use Cystoscope - Standard Deflection | | | | GA2 | BB38-01D | 30 |
| 205W11F50 | Suction Access Sheath 11/13Fr 50cm | | | | GA2 | BA35-02B | 10 |

Date: 4/17/2026 3:01 AM **Total Pallets:** **Print Sort:** Bin Location
Total Orders: 9 **Total Cartons:** **Total Pick Labels:** 1
Picker: **Pick Date:** **Batch No.:** B20260417030051
Zone: **Total Item Qty:** 263 **Pick Type:** One
Bin Category: **Conveyable:** **Schedule Ship Date:** 4/16/2026
Lift Type: **Rework:** N **Expedite:** True
Cust. Specific Attr.: **Individual:**
Consolidated: Yes **Smart Picking:** No
Pick List No.: B20260417030051One001



| Item No. | Description | Rev. | Rev. Desc. | Kit Item No. | Loc. | Bin Loc. | Qty to Pick |
|----------|----------------|------|------------|--------------|------|----------|-------------|
| AC01 | Accesory Cable | | | | GA2 | BA41-01A | 1 |